

Disclaimer: Images are for illustrative purposes.

Ultraloc Threadlocking Compound

Product Brand: Ultraloc



Size	Product Type	Product Strength	Fixture Time	Colour	Viscosity	Temperature Range	Equivalent To	sku (Part Number)
10ml Bottle	Threadlocker	Medium	4 min / 24 hrs	Blue	Medium	-50c to +150c	3440	82UTL344010
10ml Bottle	Threadlocker	Medium	10 min / 24 hrs	Blue	High	-50c to +150c	3242 242 243 7243 TM44 T43	82UTL324210
10ml Bottle	Threadlocker	Medium	10 min / 24 hrs	Blue	High	-50c to +150c	3243 242 243 7243 TM44 T43	82UTL324310
10ml Bottle	Threadlocker	Medium	20 min / 24 hrs	Green	Low	-50c to +150c	3290 290 7290 RL65 T90	82UTL329010
10ml Bottle	Threadlocker	Medium	10 min / 24 hrs	Green	Medium	-50c to +150c	3294 294	82UTL329410
10ml Bottle	Threadlocker	Low	20 min / 24 hrs	Purple	High	-50c to +150c	3222 222 7222 TM11 T22	82UTL322210
10ml Bottle	Threadlocker	High	4 min / 24 hrs	Red	High	-50c to +150c	3760	82UTL376010
10ml Bottle	Threadlocker	High	20 min / 24 hrs	Red	High	-50c to +150c	3262 262 7262 TH62 T62	82UTL326210
10ml Bottle	Threadlocker	High	10 min / 24 hrs	Red	Medium	-50c to +150c	3271 271 TM66 T71	82UTL327110
10ml Bottle	Threadlocker	High	30 min / 24 hrs	Red	High	-50c to +230c	3272 272 7272 TH68 T72	82UTL327210
10ml Bottle	Threadlocker	High	30 min / 24 hrs	Red	High	-50c to +150c	3277 277 7277 TH68 T77	82UTL327710
50ml Bottle	Threadlocker	Medium	4 min / 24 hrs	Blue	Medium	-50c to +150c	3440	82UTL344050
50ml Bottle	Threadlocker	Medium	10 min / 24 hrs	Blue	High	-50c to +150c	3242 242 243 7243 TM44 T43	82UTL324250
50ml Bottle	Threadlocker	Medium	10 min / 24 hrs	Blue	High	-50c to +150c	3243 242 243 7243 TM44 T43	82UTL324350
50ml Bottle	Threadlocker	Medium	20 min / 24 hrs	Green	Low	-50c to +150c	3290 290 7290 RL65 T90	82UTL329050
50ml Bottle	Threadlocker	Medium	10 min / 24 hrs	Green	Medium	-50c to +150c	3294 294	82UTL329450

Size	Product Type	Product Strength	Fixture Time	Colour	Viscosity	Temperature Range	Equivalent To	sku (Part Number)
50ml Bottle	Threadlocker	Low	20 min / 24 hrs	Purple	High	-50c to +150c	3222 222 7222 TM11 T22	82UTL322250
50ml Bottle	Threadlocker	High	4 min / 24 hrs	Red	High	-50c to +150c	3760	82UTL376050
50ml Bottle	Threadlocker	High	20 min / 24 hrs	Red	High	-50c to +150c	3262 262 7262 TH62 T62	82UTL326250
50ml Bottle	Threadlocker	High	10 min / 24 hrs	Red	Medium	-50c to +150c	3271 271 TM66 T71	82UTL327150
50ml Bottle	Threadlocker	High	30 min / 24 hrs	Red	High	-50c to +230c	3272 272 7272 TH68 T72	82UTL327250
50ml Bottle	Threadlocker	High	30 min / 24 hrs	Red	High	-50c to +150c	3277 277 7277 TH68 T77	82UTL327750
250ml Bottle	Threadlocker	Medium	10 min / 24 hrs	Blue	High	-50c to +150c	3242 242 243 7243 TM44 T43	82UTL3242250
250ml Bottle	Threadlocker	Medium	10 min / 24 hrs	Blue	High	-50c to +150c	3243 242 243 7243 TM44 T43	82UTL3243250
250ml Bottle	Threadlocker	Medium	20 min / 24 hrs	Green	Low	-50c to +150c	3290 290 7290 RL65 T90	82UTL3290250
250ml Bottle	Threadlocker	Medium	10 min / 24 hrs	Green	Medium	-50c to +150c	3294 294	82UTL3294250
250ml Bottle	Threadlocker	Low	20 min / 24 hrs	Purple	High	-50c to +150c	3222 222 7222 TM11 T22	82UTL3222250
250ml Bottle	Threadlocker	High	20 min / 24 hrs	Red	High	-50c to +150c	3262 262 7262 TH62 T62	82UTL3262250
250ml Bottle	Threadlocker	High	10 min / 24 hrs	Red	Medium	-50c to +150c	3271 271 TM66 T71	82UTL3271250
250ml Bottle	Threadlocker	High	30 min / 24 hrs	Red	High	-50c to +230c	3272 272 7272 TH68 T72	82UTL3272250
250ml Bottle	Threadlocker	High	30 min / 24 hrs	Red	High	-50c to +150c	3277 277 7277 TH68 T77	82UTL3277250

Disclaimer:

Dimensional data and technical information was obtained from publicly available sources and not acquired through standard agencies. It has been completed and compiled for reference purposes only; where discrepancies are found they are subject to change without notice. Bolt & Nut Australia makes no warranties or representations regarding the accuracy and validity of the compiled information and data. Contact the relevant standard authorities for accurate and detailed information.

© Bolt & Nut Australia Pty Ltd. All rights reserved.

Version: 17 Mar 2026